

Cell Scrapers & Lifters

I. APPLICATIONS

1. Cell culture
2. Cutting and scraping of surfaces

II. FEATURES

1. Materials: polyethylene
2. Gamma radiation sterilized
3. Non-pyrogenic

III. SCRAPER/LIFTER DESIGN

1. Scrapers available in 18cm and 25cm, lifter is 18cm
2. 1.8cm blade length for scraper, 2.0cm blade length for lifter
3. Pivoting blade on scraper for precise cutting and scraping
4. Individually wrapped

IV. PACKAGING

1. Bulk, 100 pieces/case or 200 pieces/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
70-1180	Cell scraper. Handle length: 18cm, blade width: 16.68	Individually wrapped, 100 pieces/case	Yes
70-1250	Cell scraper. Handle length: 240.4cm, blade width: 16.68	Individually wrapped, 100 pieces/case	Yes
70-2180	Cell lifter. Handle length: 18cm, blade width: 18.74	Individually wrapped, 100 pieces/case	Yes

VI. STERILITY

1. Products are Gamma Radiation Sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E
70-1250	Cell Scraper	240.4	16.68	9	89.0	O.D. 7.7
70-1180	Cell Scraper	180.0	16.68	9	85.0	O.D. 7.4
70-2180	Cell Lifter	180.0	18.74	7.4	-	-

